



Sentrant™

INTRODUCER SHEATH WITH HYDROPHILIC COATING




Ensure Seal

- The choice for superior hemostasis,* **Sentrant** Introducer Sheath with hydro offers optimal seal, trackability and kink resistance for enhanced procedural confidence.

Take the next step at www.medtronicendovascular.com

Innovating for life.



➤ Compatible with Reliant Balloon and Archer Guidewire

PRODUCT ORDERING INFORMATION

Model Number	Sheath Size (F)	Working Length (cm)
SENSH1228W	12	28
SENSH1428W	14	28
SENSH1628W	16	28
SENSH1828W	18	28
SENSH2028W	20	28
SENSH2228W	22	28
SENSH2428W	24	28
SENSH2628W	26	28

Superior Hemostasis*

- Optimal seal versus competitor sheaths


Enhanced Trackability

- Hydrophilic coating and flexibility for easier tracking through tortuous calcified iliacs

Procedural Confidence


- Reinforced coil for kink resistance
- Radiopaque marker band for easy visibility
- Dilator locking mechanism secures desired position

SUPERIOR LEAK RESISTANCE VERSUS COOK CHECK-FLO PERFORMER AND GORE DRYSEAL†



* Bench Test Data on file at Medtronic, Inc. Test data not indicative of clinical performance. Bench Test compared Cook Check-Flo Performer and Gore DrySeal 12F and 18F to Sentrant 12F and 18F.
 † Medtronic sheath had zero leakage on 0.035" guidewire.
 ‡ Both Sentrant 12F and 18F were tested and demonstrated superior kink resistance. 18F results shown above.

SUPERIOR KINK RESISTANCE VERSUS COOK CHECK-FLO PERFORMER AND GORE DRYSEAL‡



Sentrant™ Introducer Sheath with Hydrophilic Coating

Important Information: Prior to use, refer to the *Instructions for Use* supplied with these devices for indications, contraindications, suggested procedure, warnings and precautions.
Indications for Use: The Medtronic Sentrant Introducer Sheaths with Hydrophilic Coating are intended to provide a conduit for the insertion of diagnostic or endovascular devices into the vasculature and to minimize blood loss associated with such insertions.
Potential Adverse Events: The following adverse events or complications may be associated with the use of the Sentrant introducer sheath: allergic response to materials; blood loss, bleeding, or hematoma; embolization (micro or macro) with transient or permanent ischemia; infection; vascular trauma (e.g., dissection, rupture, perforation, or tear); and death.
Warnings: This device should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device.
CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a licensed healthcare practitioner. See package inserts for full product information.

➤ www.medtronic.com
www.medtronicendovascular.com

Medtronic Vascular, Inc. 3576 Unocal Place Santa Rosa, CA 95403 USA	Product Services Support Center Tel: 888.283.7868 Fax: 800.838.3103	CardioVascular LifeLine Customer Support Tel: 877.526.7890 Tel: 763.526.7890
---	--	---

